Chapter 3 – Week 1 Parts A,B,C

Purification of Lactate Dehydrogenase (LDH)

Purpose:

1. Learn basic techniques for protein purification

2. Prepare crude extract

3. Run enzyme assay

Protein Purification Process


Tis	sues or Cells	 Disruption by: Ultrasonic vibration French press Blending 		
Hom	ogenate or Cell Extract	 Separation by centrifugation Pellet – Cell Debris Supernatant – Cell Extract 		
Imp	oure Proteins	 Separation by column chromatography or ionic strength Ammonium Sulfate Precipitation Affinity Chromatography 		
	Purified	Separation by column chromatographyGel-Filtration ChromatographyIon-Exchange Chromatography		
Pure Protein				

Centrifugation

- Uses centrifugal force applied to sample to separate complex mixtures
 - Breaking down of cells into various components
- Largest units experience largest centrifugal force, moving rapidly to bottom of tube
- Speed is observed in RPM revolutions per minute
- Can compare between different centrifuges using RCF – relative centrifugal force:

 $\begin{aligned} \text{RCF} &= [r(2\pi N)^2]/g = (1.118 \times 10^{-5})(r_{\text{cm}})(N_{\text{RPM}})^2 \\ g &= \text{gravitational acceleration} = 9.8 \text{ m/s}^2 \\ r &= \text{radius of rotor} \\ N &= \text{rotational speed in RPM} \end{aligned}$






Lactate Dehydrogenase (LDH)

LDH catalyzes the last step of anaerobic glycolysis

 ⁰
 ¹
 ⁰
 ¹
 ⁰
 ¹
 ¹




- Known kinetic parameters and crystal structure
- Multiple forms of LDH found in different tissues

 Isozymes
 - Each isozyme has slightly different kinetic and structural properties, but same function and overall structure

Cofactors of LDH


- LDH uses pyridine nucleotide coenzymes as cofactors to transfer reducing equivalents
 - H⁺ and a pair of electrons
- NAD⁺ Nicotinamide adenine dinucleotide \rightarrow Oxidized form
- NADH \rightarrow Reduced form




- Pyridine nucleotides have characteristic absorbance spectra
 - NAD⁺: No absorbance at 340 nm
 - NADH: $\lambda_{max} = 340 \text{ nm}, \epsilon = 6210 \text{ M}^{-1} \text{ cm}^{-1}$

How to Measure if LDH is Present?

Enzymes are catalysts – speed up the rate of the reactions without being consumed



- We can measure:
 - Rate of consumption of reactants (Pyruvate or NADH)
 - Rate of formation of products (L-Lactate or NAD⁺)
- NADH has a visible absorbance at 340 nm can follow its consumption during reaction
- Always want to monitor the initial rate because concentrations change during the reaction



- Purify Beef Heart or Skeletal Muscle
- Absorption Spectra & Calculate Extinction Coefficient
- Activity Assay
- Protein Concentration Dye Binding Assay

Keep everything on ice! Especially LDH extract!

- Purify Beef Heart or Skeletal Muscle
 - Blend with buffer, measure volume of homogenate
 - Centrifuge 15,000 rpm for 15 minutes
 - Discard 1P fraction, measure total volume of 1S fraction
 - Save ~1 mL of 1S fraction in separate aliquot for activity assays and dye-binding
 - Save rest of 1S fraction for Week 2

Record everything! Masses, volumes, concentrations etc! When in doubt, don't throw it out!

- Absorption Spectra & Calculate Extinction Coefficient
 - Obtain spectra of NAD⁺ and NADH with UV-Vis Spectrometer
 - Calculate the apparent extinction coefficient for your bench top spectrophotometers
 - Given extinction coefficients are for square cuvettes with 1 cm path length
 - Measure A₃₄₀ of NADH with your spectrophotometer and UV-Vis spectrophotometer and calculate:

 $\varepsilon_{\lambda}(apparent) = [\varepsilon_{\lambda}(UV)][A_{\lambda}(test tube)/A_{\lambda}(UV cuvette)]$

- Activity Assay
 - Monitors $A_{\rm 340\,nm}$ to observe of the conversion of NADH to NAD+ as the reaction proceeds
 - Requires cocktail with all other reagents except enzyme and water
 - Make enough for 10 assays (See Recipe, p. 73)

	Reagent	Volume
"The	1 M KPO ₄ , pH 7.4	200 μl
	6 mM Sodium pyruvate	300 μl
	1 mM NADH	300 μl
	DI Water	2.15 ml
	Diluted Enzyme Solution	50 μl
	Total Volume	3.00 ml

- Will need to dilute enzyme substantially before starting assays
 - Typically 1:200 or 1:400 dilutions of enzyme are required

- Activity Assay
 - To run assays:
 - Use LDH Kinetics program
 - Blank spectrometer
 - Add cocktail and water to cuvette
 - Add enzyme, mix, and start run
 - Want to look at linear portion of graph
 - Readout is a RATE: $-\Delta A_{340}/min$
 - Why Negative?

•Need at least 4 dilutions in range ΔA_{340} /min -0.05-0.25 •Need table of time vs. A_{340} for 1 run

- Protein Concentration Dye Binding Assay
 - Do only if time allows, otherwise save for week 2
 - You can use your standard curve from Lab 1 if it goes through the origin and is a good fit
 - When in doubt, make a new standard curve
 - Record protein concentration for 1S fraction
 - Dilute your protein with water!
 - Blank?
 - Use protein concentration with activity for calculations

Activity Calculation

Activity = Units = μ mol of Substrate Consumed or Product Formed / min Activity concentration (ΔC) = [Activity] = $\mu mol/min^*ml$ = Units/ml From Beer's Law: $A = (\varepsilon)(I)(C)$ We are looking at Kinetic Rates: $\Delta A = (\epsilon)(I)(\Delta C)$ Apparent extinction coefficient, ε_{app} , in mM⁻¹ combines ε and I $\Delta A_{340/\text{min}} = (\epsilon_{app} \text{ in mM}^{-1}) (\Delta C)$ [Activity] = $\Delta C = \Delta A_{340/min} / \epsilon_{app}$ in mM⁻¹ $[Activity] = (0.05/min)/(6.21 \text{ mM}^{-1}) = 0.0081 \text{ mM}/min =$ 0.0081 µmol/min*ml = 0.0081 units/ml $mM = \mu mol/ml$ Enzyme Activity in the Assay $mM^{-1} = ml/\mu mol$

Activity Calculation

[Activity] = Units/ml = µmol of Substrate Consumed or Product Formed

min * ml

You must account for the dilutions of your protein!

 $[Activity_{Undiluted}] = (\Delta C)(Total Volume of Assay)(Dilution Factor)$ (Volume of enzyme used in assay) $[Activity_{Undiluted}] = (0.0081 \text{ units/ml})(3.0 \text{ ml})(400) = 193 \text{ units/mL}$ (0.05 ml)

More Enzyme Calculations

- Total Activity = (Activity)(Total Volume) = Units/ml* ml = Units
- Protein = Mass Protein/Volume Extract = mg/ml
- Total Protein = (Protein)(Total Volume) = mg/ml* ml = mg
- Specific Activity = Total Activity/Total Protein = Units/mg
- % Yield = Total Activity in Given Step
 Total Activity in Crude Extract x 100

Remember to account for the dilutions of your protein!

For needed calculations, see purification table, p. 86