

# Announcements

- Quiz on Chapter 5a today, Chapter 5b and 5c quiz next week
- Rest of semester schedule:

	Discussion Dates	Lab Dates	Lab Due Dates
Chapter 5a		Mon. 11/5	
Chapter 5b	11/5 – 11/7	11/7 – 11/12	
Chapter 5c	11/12 – 11/14	11/14 – 11/19	11/28 – 12/3 in lab
Chapter 6a	11/26 – 11/28	11/28 – 12/3	
Chapter 6b	12/3 – 12/5	12/5 – 12/10	12/12, All Sections!
Lab Exam	12/10 – 12/12		

# Chapter 5: Structural Characterization of LDH

## Purpose:

Learn how to:


- 1) Pour an SDS Polyacrylamide gel
- 2) Prepare samples for SDS-PAGE
- 3) Determine purity and subunit MW of your LDH

# Chapter 5: Overview

- **Week 1: Gel Filtration Chromatography**
  - Separate LDH from standard proteins by native size
  - Determine Native MW
- **Week 2: Denaturing Polyacrylamide Gel Electrophoresis (SDS-PAGE)**
  - Separate LDH from other proteins by subunit size
  - Determine Subunit MW
- **Week 3: Native Electrophoresis (Zymograms)**
  - Confirm the quaternary structure of LDH from Weeks 1 & 2

# Separating Proteins by Electrophoresis

- **Applied electric field separates proteins**
- **Many factors effect rate:**
  - Isoelectric points
  - Titration curves
  - Molecular weight
  - Hydrodynamic properties




- **Native Electrophoresis**
  - Separation of full protein by size and charge
  - Maintains tertiary and quaternary structure
- **Denaturing Electrophoresis**
  - Denature of protein into subunits
  - Separation of those denatured subunits by size
- **2D Electrophoresis**
  - Separation of protein by pI in first direction
  - Separation by size in second direction

# Polyacrylamide Gels


- Matrix for separation
- Acts as homogeneous support to prevent diffusion of proteins out of gel
- Acts as **molecular sieve** slowing the migration of proteins *in proportion to their charge-to-mass ratio*
- Concentration can be varied to separate wide range of proteins
- Speed of separation is related to **electrophoretic mobility**

# Polyacrylamide Gel Formation

Acrylamide (monomer)




+



N,N'-Methylene-Bis-acrylamine (crosslinker)

Free Radical  
Chain Reaction

Initiated by TEMED  
and Persulfate



Need to remove any un-polymerized monomer – it can react with protein functional groups!


Polyacrylamide Gel Composition

# Denaturing Polyacrylamide Gel Electrophoresis – SDS-PAGE

- Method to separate proteins by molecular weight of denatured subunits
- Treat protein with strong denaturant (SDS) and sulfhydryl compound ( $\beta$ ME or DTT)
  - Unfolds polypeptide chains into random coils
  - Coats protein surface to give uniform negative charge
  - Prevents disulfide bonds from forming between subunits
- Denatured subunits can be permanently separated with SDS containing buffers
  - Large excess of SDS
  - SDS binds to any protein at 1.4 g SDS/g of protein
    - Gives all proteins same charge/mass ratio

# SDS-PAGE

- Discontinuous (Disc) electrophoresis
  - Layers of gel for better separation based on pH changes (p. 165-8)
- **Running Gel** – bottom layer, high pH (8.8), higher % acrylamide
- **Stacking Gel** – top layer, lower pH (6.8), lower % acrylamide, where the wells are located
- Stacking layer for loading samples & stacking condition, running layer for separation
- Different acrylamide percentages in gel affect separation



- **Higher % gel** – better separation of small proteins
- **Lower % gel** – better separation of large proteins





# Relative Mobility

## Linear relationship between mobility of a protein and Log MW

- Measure distance from top of gel to solvent front
- Measure distance from top of gel to protein band interest

$$R_m = \frac{\text{Distance to protein band on interest}}{\text{Distance to Solvent Front}}$$



**Protein Standards  
used in Lab 5 on p.  
156 of Lab Manual**

# Procedure: Chapter 5-Week 2

- Prepare Samples
- Assemble Gel Apparatus
- Load Samples
- Run Gel
- Pour Gels (for the next lab section)
- Staining and Visualization

# Procedure: Chapter 5-Week 2

## • Prepare Samples

- Crude Extract – 30-50  $\mu\text{g}$
- 3P-Dialyzed – 15-30  $\mu\text{g}$
- Purified LDH – 2-5  $\mu\text{g}$
- SDS-PAGE Standards –  
5  $\mu\text{l}$  aliquot, ready to load, 1 per gel

## • Total Sample Volume = 20 $\mu\text{l}$

- 2  $\mu\text{l}$  10X Sample Buffer (thick blue liquid)
- Up to 18  $\mu\text{l}$  of sample to put at the appropriate concentration of protein
- Water to bring total volume to 20  $\mu\text{l}$

## • Denature samples 1-2 min at $\sim 90^\circ\text{C}$

- ❖ Use concentration from Dye-Binding from Chapter 3!
- ❖ Calculate volumes before you come to the lab!

### **Standards:**

Bio-Rad Broad Range  
See Table p. 156  
Marked with ¶ Symbol

# Procedure: Chapter 5-Week 2

- **Assemble Gel Apparatus:**

- 2 groups/gel, 2 gels/gel box
- Pictures p. 135-136 and demo!

- **Load Samples:**

- Use gel loading tips
- Yellow loading guides can help line up wells
- Write down your loading scheme!

- **Run Gel:**

- 50 V until through the stacking gel, 200V after that
- Run until blue tracking dye reaches bottom of gel
- Total time ~ 1 hr

# Procedure: Chapter 5-Week 2

- Pour Gels:

- Put together short plate and spacer plate (1.5 mm)
- Lock together with casting clamp and stand (pictures p. 133)
- Check plates for leaks with water
- Prepare Running Gel – Recipe p. 134 **(100  $\mu$ l APS)**
  - Pour or pipet into plate, save extra and see when it polymerizes
  - Overlay with isobutanol
- Prepare Stacking Gel – Recipe p. 134 **(40  $\mu$ l APS)**
  - Pour or pipet into plate, save extra and see when it polymerizes
  - Insert comb to form wells
- Put poured gel in buffer at 4°C for the next section

# Procedure: Chapter 5-Week 2

- Staining and Visualization:

- Remove gel from plates
- Nick a corner so you know the orientation
- Stain with Coomassie Brilliant Blue R-250
- Destain with 1:5:5 Acetic Acid/Methanol/Water solution
- Image on gel doc
- Calculate  $R_m$  values for standards and unknown protein bands
- Make graph of Log MW vs  $R_m$



# Calculate $R_m$ Values

## Make Plot of Log MW vs. $R_m$

### Linear relationship between mobility of a protein and Log MW

- Measure distance from top of gel to solvent front
- Measure distance from top of gel to protein band interest

$$R_m = \frac{\text{Distance to protein band on interest}}{\text{Distance to Solvent Front}}$$



**Protein Standards  
used in Lab 5 on p.  
156 of Lab Manual**