

Understanding secondary antibodies

Fragment antigen binding antibodies and isotypes





www.abcam.com/secondary_antibody

Antibody structure and F(ab) antibodies

The light chain (LH) folds into a variable domain (VL) and a constant domain (CL) whereas the heavy chain is composed of one variable domain (VH) and three (IgG and IgA) or four constant domains (IgE).





The **F(ab)** fragment is an antibody structure that still binds to antigens but is monovalent with no Fc portion. An antibody digested by the enzyme papain yields two F(ab) fragments of about 50 kDa each and an Fc fragment. In contrast, **F(ab')₂ fragment** antibodies are generated by pepsin digestion of whole IgG antibodies to remove most of the Fc region while leaving intact some of the hinge region. F(ab')₂ fragments have two antigen-binding F(ab) portions linked together by disulfide bonds, and therefore are divalent with a molecular weight of about 110 kDa.

Secondary antibodies

F(ab) fragments

Host	Target	Conjugate	Applications	www.abcam.com/
Donkey	Goat IgG - H&L	HRP	Dot, ELISA, EM, IHC, WB	ab6667
	Goat IgG - H&L	TRITC	Flow Cyt, IM	ab6522
	Goat IgG - H&L	TR	Flow Cyt, IM	ab6523
	Goat IgG - H&L	Biotin	Dot, ELISA, IHC, IM, WB	ab6578
Goat	Human InG - H&I	1 4nm Gold	EM WB	ab30826
Cour	Mouse IgG - E(ab') pre-adsorbed		Flow Cyt. ICC/IE IHC-P	ab51379
	Mouse IgG - H&L	Intro	110w Oyt, 100/11, 1110-1	ab6668
	Mouse IgG - H&L	HRP	Dot, ELISA, EM, IHC, WB	ab6823
	Mouse IgG - H&L	TRITC	Flow Cyt, IM	ab6670
	Mouse IgG - H&L	FITC	Flow Cyt, IM	ab6669
	Rabbit IgG - H&L			ab6824
	Rabbit IgG - H&L	Biotin	Dot, ELISA, IHC, IM, WB	ab7055
	Rabbit IgG - H&L	HRP	Dot, ELISA, EM, IHC, WB	ab7171
	Rat IgG - H&L	Biotin	Dot, ELISA, IHC, IM, WB	ab7176

F(ab')₂ fragments

Host	Target	Conjugate	Applications	www.abcam.com/
Donkey	Goat IgG - H&L, pre-adsorbed Rabbit IgG - H&L, pre-adsorbed	PE PE	Flow Cyt, IHC-Fr, IHC-P, IM Flow Cyt, ICC/IF, IHC-P, IM	ab7004 ab7007
Goat	Mouse IgG - (Fab)' ₂ , pre-adsorbed Mouse IgG - Fc, pre-adsorbed Mouse IgG - H&L, pre-adsorbed Mouse IgG+IgM+IgA - H&L, pre-adsorbed Mouse IgG+IgM+IgA - H&L, pre-adsorbed Mouse IgM - mu chain Mouse IgM - mu chain Mouse IgM - mu chain Mouse IgM - mu chain Rabbit IgG - (Fab)' ₂ , pre-adsorbed Rabbit IgG - (Fab)' ₂ , pre-adsorbed	TR PE PE HRP FITC Biotin HRP AP FITC HRP TRICT	Flow Cyt, IM Flow Cyt, IM Flow Cyt, IHC-Fr, IHC-P, IM Dot, ELISA, EM, IHC, WB Flow Cyt, IF Dot, ELISA, Flow Cyt, IHC, IM Dot, ELISA, IHC, WB Flow Cyt, IM Dot, ELISA, EM, IHC, WB Flow Cyt, IM	ab5884 ab5881 ab7002 ab6006 ab5999 WB ab5929 WB ab5930 ab5931 ab5926 ab6112 ab6110
Rabbit	Goat IgG - H&L Goat IgG - H&L	HRP FITC	Dot, ELISA, EM, IHC, WB Flow Cyt, ICC/IF	ab5755 ab47846
Sheep	Mouse IgG - H&L , pre-adsorbed	Cy3®	IHC-P	ab50502

See more species and isotype reactivity at www.abcam.com/search_secondary

Engineered antibodies



- A single-chain variable fragment (scFv) is a powerful tool in phage display, as it can express the antigen-binding domain as a single peptide. ScFv is also used in FACS, IHC and as antigen-binding domains of artificial T cell receptors.
- **Diabody** plays a role in the recruitment of antibody effector functions and cytotoxic T-cell responses.
- A **single domain antibody** (sdAb) consists of a single monomeric variable antibody domain and is involved in pharmaceutical applications.
- A **trifunctional antibody** is a monoclonal antibody with binding sites for two different antigens, typically CD3 and a tumor antigen. Trifunctional antibodies are considered as a type of bispecific monoclonal antibody.
- **Bi-specific T-cell engagers** (BiTE[®]) is a class of bispecific monoclonal antibody that direct a host's T cells cytotoxic activity against diseased cells.

F(ab) and Fc Receptors



names	FcγRICD64	FcyRIIB2CD32	FcɛRllCD23
eptor I	FcyRIIACD32	FcyRIIIACD16a	FcαRI CD89
Fc rec	FcyRIIB1CD32	FcyRIIIBCD16b	

F(ab) and $F(ab')_2$ antibodies eliminate non-specific binding between Fc portions of antibodies and Fc receptors on cells (such as macrophages, dendritic cells, neutrophils, NK cells, B cells...)

Ig structural differences



Isotype: Distinct forms of light and heavy chains which are present in all members of a species. Kappa and lambda are isotypes of light chains; mu, delta, gamma, alpha and epsilon are isotypes of heavy chains.



Allotype: Allelic variants within the constant region of the immunoglobulin light or heavy chains. Of a given isotype, members of a species differ in function to the particular alleles they have received from their parents. **Idiotype:** Antigenic specificity of a particular monoclonal immunoglobulin.

Secondary antibodies

Anti-IgG H&L

DonkeyGoatHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6885GoatMouseFITCFlow Cyt, ICC/IF, IHC-Fr, IHC-P, IM, WBab6785MouseHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6786RabbitAPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6722RabbitBiotinDot, ELISA, IHC-Fr, IHC-P, IM, WBab6720RabbitBiotinDot, ELISA, IHC-Fr, IHC-P, IM, WBab6720RabbitCy3*Flow Cyt, ICC/IF, IHC-Fr, IHC-Pab6639RabbitCy5*Flow Cyt, ICC/IF, IHC-Fr, IHC-Pab66741RabbitFITCFlow Cyt, ICC/IF, IHC-Fr, IHC-P, IMab6717RabbitHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6721RabbitRabbitHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6721RabbitGoatHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6726MouseTRFlow Cyt, ICC-Fr, IHC-P, IM, WBab6726RatHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6726RatHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6734SheepHRPDot, ELISA, IHC-Fr, IHC-P, IM, WBab6734

Anti-IgM

Host	Target	Applications	Datasheet www.abcam.com/
Goat	Human Human Mouse Mouse (Biotin) Rabbit	Dot, ELISA ELISA, ID, Ie, WB ELISA, ID, Ie, WB ELISA ELISA ELISA, ICC, WB	ab8496 ab9160 ab9167 ab49760 ab2754
Mouse	Sheep (clone [G9])	Conjugation, ELISA	ab19991
Rabbit	Human Mouse	Dot, ELISA ELISA, Ie, WB	ab8505 ab9175
Rat	Rabbit (clone [SB87a])	ELISA, Flow Cyt	ab80301

Anti-IgA

Host	Target	Applications	Datasheet www.abcam.com/
Chicken	Human	ELISA, ID, Ie, WB	ab9151
Goat	Human	Dot, ELISA	ab8499
	Human	ELISA, ID, Ie, WB	ab9157
	Mouse	ELISA, ID, Ie, WB	ab9161
	Rabbit	ELISA, ICC, WB	ab2758
Mouse	Human (clone [hlgA1])	ELISA, RIA	ab433
Rabbit	Human	Dot, ELISA	ab8508
	Mouse	ELISA, ID, Ie, WB	ab9170

Enzymatic detection methods

Enzyme	Substrates	Applications	Advantages	Disadvantages
Horseradish Peroxidase (HRP)	Chromogenic, soluble (TMB, ABTS, OPD)	ELISA	Easy to use	Light sensitive coloration
	Chromogenic, precipitating (CN, AEC, DAB,)	WB, SB, IHC	Easy to use	Background in blood samples and some other tissues Staining stability lower than AP
	Fluorogenic (ADHP/resorufin) Chemiluminescent	ELISA	High sensitivity	Need fluorescence equipment
	Luminol	WB, SB, IHC	High sensitivity	Need radiographic equipment or light scanner
Alkaline Phosphatase (AP)	Chromogenic, soluble (pNPP)	ELISA	Linear kinetic Often more sensitive than HRP	Unstable
	Chromogenic, precipitating (BCIP/NBT,)	WB, SB, IHC	Staining stability higher than HRP	Interference with nuclear counterstain
	Fluorogenic (4-MUP)	ELISA, IHC	Sensitivity	Need fluorescence equipment

WB: Western Blot, SB: Southern Blot, IHC: Immunohistochemistry, ELISA: Enzyme Linked Immunosorbent Assay

How to choose the right seconda

Guide to help you choose the most appropriate secondary antibody for your application.

1. What is the host species of the primary antibody?

The secondary antibody is directed against the species of the primary antibody. If you use a primary antibody raised in rabbit, you will need an antirabbit secondary antibody raised in a species other than rabbit.

2. What do I need to know about the isotype of the primary antibody?

The secondary antibody has to be directed against the isotype of the primary antibody. Polyclonal primary antibodies are generally raised in rabbit, goat, sheep or donkey and are an IgG isotype. The secondary antibody will typically be an anti-IgG H&L antibody.

Monoclonal primary antibodies are commonly raised in mouse, rabbit and rat. For example, if the primary monoclonal antibody is a mouse IgG1, you will need an anti-mouse IgG or a less specific F(ab) fragment anti-mouse IgG.

Human immunoglobulin classes, subclasses, types and subtypes:

- Classes or isotypes: IgG (γ heavy chains), IgM (μ), IgA (α), IgE (ϵ), IgD (δ)
- Subclasses: IgG1 (γ1 heavy chains), IgG2 (γ2), IgG3 (γ3), IgG4 (γ4), IgA1 (α1), IgA2 (α2)
- Types: κ light chain, λ light chain
- Subtypes: λ1, λ2, λ3, λ4

Other type of reactivities:

- · Polyvalent antibodies react with all classes
- Anti-Fc or heavy chain (α , δ , ϵ , γ , and μ) antibodies react with heavy chain only
- Anti-F(ab) or whole molecule antibodies react with heavy and light chains independently of the class
- Anti-light chain (κ and λ) antibodies react with all classes since all classes use the same κ and λ light chains

3. Do I need an enzymatic or fluorescent detection?

The type of conjugation is application dependent.

For enzymatic and biotin detection, e.g. in WB or ELISA, we suggest a secondary antibody conjugated to HRP, AP or biotin. Both avidin and streptavidin bind very strongly to biotin and enable signal amplification,

ry antibody?

regardless of the host species of the antibody.

If a laser light is involved, e.g. in Flow Cytometry, ICC/IF or IHC, we suggest fluorescent detection with a secondary antibody conjugated to a fluorochrome.

4. Do I need a pre-adsorbed secondary antibody?

We usually recommend using a secondary antibody, pre-adsorbed with serum, for western blotting, of immunoglobulin-rich tissues and cells. Pre-adsorbed secondary antibodies are less likely to interact with endogenous immunolgobulins and consequently may reduce non-specific background. The secondary antibody should be pre-adsorbed against the same species as the sample on which the detection is performed. For example, a human pre-adsorbed antibody will be required for detection in human tissue.

5. Do I need an affinity purified antibody or IgG fraction?

The advantage of using affinity purified antibodies or IgG fractions will depend on the type of binding expected. Affinity purified antibodies give the lowest amount of non-specific binding whereas IgG fractions contain high affinity antibodies. Indeed, during an affinity purification, high affinity antibodies stay fixed on the matrix and cannot be eluted.

6. Is it necessary to use a F(ab) or F(ab')₂ fragment antibody?

F(ab) and $F(ab')_2$ fragment antibodies eliminate non-specific binding between Fc portions of antibodies and Fc receptors on cells (such as macrophages, dendritic cells, neutrophils, NK cells and B cells) and penetrate tissues more efficiently due to their smaller size. As fragment antibodies do not have Fc portions, they do not interfere with anti-Fc mediated antibody detection.

7. Do I need an anti-IgG H&L, anti-Iight chain or anti-F(ab')₂ secondary antibody?

Our secondary antibodies are supplied in different formats:

- Anti-IgG H&L antibodies react with both heavy and light chains of IgG subclasses
- Anti-light chain antibodies react with the light chain of primary antibodies which is the same among all classes
- Anti-F(ab')₂ secondary antibodies react with the F(ab')₂ portion of the primary antibody

For more information and links to all our downloads please visit: www.abcam.com/secondary_antibody

Fluorescent detection methods

Spectral properties of DyLight® fluorochromes

Fluorochrome Color		Ex/Em	ε (Μ¹c	m ⁻¹) Spectrally equivalent dyes
DyLight [®] 488		493/518	70K	Alexa Fluor 488°, Fluorescein, FITC, Cy2°
DyLight [®] 549		562/576	150K	Alexa Fluor 546° or 555°, Cy3°, TRITC
DyLight [®] 594		593/618	80K	Alexa Fluor 594°, Texas Red [®]
DyLight [®] 649		654/673	250K	Alexa Fluor 647°, Cy5°

Ex/Em: Excitation and emission wavelength in nanometers (+/- 4nm) ϵ : Molar extinction coefficient at the absorption maximum



More information at www.abcam.com/DyLight

DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Abpromise[®] Our promise to you

100% technical support

- Thinking of purchasing a secondary antibody but have some technical questions?
- · Unsure of buffers or incubation times?
- Having difficulty with a product? ... Our scientific experts are here to help.

All Abcam secondary antibodies are guaranteed to work as specified on the datasheet.

- 100% Scientific and Customer Support for any product you buy from Abcam or one of its authorized distributors
- We guarantee our products work in the tested species and applications as stated on the datasheet
- We will replace or refund products not performing as stated on the datasheet if reported within 120 days of purchase



Copyright © 2010 Abcam, All Rights Reserved. The Abcam logo is a registered trademark. All information/detail is correct at time of going to print.