### **Announcements**

• Chapter 3 and 11 are both due at the end of Chapter 4 lab


# **Chapter 4: Enzyme Kinetics**

#### **Purpose:**

- A) Re-assay LDH activity using Chapter 3 cocktail (\*containing potassium phosphate buffer)
- B) Assay LDH activity using varying concentrations of one substrate (pyruvate)  $\rightarrow$  \*Tris buffer cocktail
  - To determine  $K_m \& V_{max}$  of pyruvate
- C) Characterize the effect of an unknown inhibitor on enzyme activity  $\rightarrow$  \*Tris buffer cocktail
  - Solve for  $K_1$


# Review: LDH catalyzes the last step of anaerobic glycolysis

- Multiple forms of LDH found in different tissues
  - → Isozymes
  - Each isozyme has slightly different kinetic and structural properties, but same function and overall structure




#### **Michaelis-Menten Kinetics**

$$E + S \xrightarrow{k_1} ES \xrightarrow{k_2} E + P$$




#### Lineweaver-Burk Manipulation

#### **Double-Reciprocal Plot**




- Competitive Inhibition Binds to free enzyme only & competes with substrate for active site
- Uncompetitive Inhibition Binds to distinct site from substrate active site and binds only to ES complex
- Pure Non-Competitive Inhibition Binds to a distinct site on the free enzyme or ES complex that decreases overall activity  $(K_1 = K_1')$
- Non-Competitive Inhibition (Mixed) Binds to free enzyme or ES complex  $(K_1 > K_1')$  or  $K_1 < K_1'$ )



- Competitive Inhibition Binds to free enzyme only & competes with substrate for active site
- Uncompetitive Inhibition Binds to distinct site from substrate active site and binds only to ES complex
- Pure Non-Competitive Inhibition Binds to a distinct site on the free enzyme or ES complex that decreases overall activity  $(K_1 = K_1')$
- Non-Competitive Inhibition (Mixed) Binds to free enzyme or ES complex  $(K_1 > K_1')$  or  $K_1 < K_1'$ )


- Competitive Inhibition Binds to free enzyme only & competes with substrate for active site
- Uncompetitive Inhibition Binds to distinct site from substrate active site and binds only to ES complex
- Pure Non-Competitive Inhibition Binds to a distinct site on the free enzyme or ES complex that decreases overall activity  $(K_1 = K_1')$
- Non-Competitive Inhibition (Mixed) Binds to free enzyme or ES complex  $(K_1 > K_1')$  or  $K_1 < K_1'$ )

- Competitive Inhibition Binds to free enzyme only & competes with substrate for active site
- Uncompetitive Inhibition Binds to distinct site from substrate active site and binds only to ES complex
- Pure Non-Competitive Inhibition Binds to a distinct site on the free enzyme or ES complex that decreases overall activity  $(K_1 = K_1')$
- Non-Competitive Inhibition (Mixed) Binds to free enzyme or ES complex  $(K_1 > K_1')$  or  $K_1 < K_1'$ )




- Competitive Inhibition Binds to free enzyme only & competes with substrate for active site
- Uncompetitive Inhibition Binds to distinct site from substrate active site and binds only to ES complex
- Pure Non-Competitive Inhibition Binds to a distinct site on the free enzyme or ES complex that decreases overall activity  $(K_1 = K_1')$
- Non-Competitive Inhibition (Mixed) Binds to free enzyme or ES complex  $(K_1 > K_1')$  or  $K_1 < K_1'$


# **Competitive Inhibition**




# **Uncompetitive Inhibition**




# Pure Non-Competitive Inhibition



# Non-Competitive Inhibition (Mixed)





See pp. 99-100 for equations

### **Chapter 4A-B: Procedure**

| Reagent   | Chapter 3 recipe    | Chapter 4 recipe |
|-----------|---------------------|------------------|
| Enzyme    | LDH                 | LDH              |
| Cofactor  | NADH                | NADH             |
| Substrate | Pyruvate            | Pyruvate         |
| Buffer    | Potassium phosphate | Tris             |
| Diluent   | dH <sub>2</sub> O   | $dH_2O$          |

- Re-assay LDH using Chapter 3 recipe and write down activity concentration
- Make new cocktail with Tris-Buffer pH 8.2 Cocktail A
  - Phosphate acts as a mild inhibitor on LDH
- Perform activity assays where you vary [pyruvate] without inhibitor
  - Starting  $\Delta A_{340}$ /min = 0.02-0.04 for lowest [pyruvate] (\*correction)
  - Dilute appropriately to get in range

### **Chapter 4C: Procedure**

- Make new cocktail with Tris-Buffer pH 8.2 and inhibitor – Cocktail B
  - Make sure to write down letter and concentration of your assigned inhibitor
- Perform activity assays where you vary [pyruvate] in presence of the inhibitor
  - Rates with inhibitor < Rates of uninhibited reactions</li>

Make sure to prepare data tables p. 106-7
BEFORE LAB!

Include all cocktail recipes in your notebook!

### **Chapter 4: Lab manual typos**


- Cocktail tables (pg 104 & 105) --
  - For Na<sup>+</sup> pyruvate, volume range should be 0.025 0.25 mL
  - For H<sub>2</sub>O, volume range should be 0.475 0.25 mL
- Part B.1 & Table 1 in notebook section
  - For your initial trials, the lowest pyruvate concentration should be 0.05 mM, not 0.1 mM
- ΔA/min range for initial trials
  - For your initial trials, use a ΔA/min range between 0.02 and 0.04 at the lowest pyruvate concentration

# Lab Notebook: Chapter 4

- Raw Data for uninhibited and inhibited LDH
- Calculation of rates in mM/min:

$$\left(\frac{\left(\frac{\Delta A_{340}}{min}\right)}{\left(\varepsilon_{app\ in\ mM^{-1}}\right)}\right)\left(\frac{(3\ mL\ total\ volume)(Dilution\ Factor)}{(0.1\ mL\ enzyme\ used)}\right) = Rates\ in\ mM/min$$

 Michaelis-Menten and Lineweaver-Burk Plots for uninhibited and inhibited LDH




- Calculation of K<sub>M</sub> and V<sub>max</sub> Show calculations!
- Calculation of K<sub>1</sub> for <u>your</u> type of inhibition

# Lab Notebook: Chapter 4

• To calculate  $K_{l}$  you will need your  $V_{max}$ ,  $K_{m}$  and inhibitor concentration values, depending on your type of inhibitor

Plot your data first to figure out which type of inhibitor you have



### **Chapter 4:**

#### Before the lab period, you should have:

- ✓ Completed your prelab
  - ✓ Title, date, introduction, procedures
  - ✓ Be sure to account for all corrections and changes in lab manual

#### At the end of lab, you should have:

- ✓ Re-assayed your LDH using Chapter 3 cocktail
- ✓ Collected LDH data for varying pyruvate concentrations
- ✓ Recorded your inhibitor letter and concentration
- ✓ Collected LDH data for varying [pyruvate] with inhibitor present

# Questions?