



## Chapter 10BC: Lipids and Membranes Note that you need to do appropriate back-calculations and account for **Objectives** theoretical masses, that is, "what if you had not washed any portion of the membrane?" for your unwashed fraction and "what if you had processed all of your membrane for washing?" for your supernatant fraction and washed fractions. Think carefully when doing your post-lab write up! From your RBC membrane preps you will determine: Weight Data Fractions Weight of Weight of Weight of Weight of protein\*\*\* cholesterol\*\*\*\* fractions\* lipid\*\* Unwashed Washed SF \* Equal to sum of (weights of lyophilized preps (week 2) - weight of tared tubes) (week 1) \*\* Equal to weight of dried extracted lipids (week 2) - weight of tared tube (week 2) \*\*\* Equal to conc. of protein in aliquots from standard curve (week 2) x volume of each fraction (week 1) x 1/fraction of prep aliquoted \*\*\*\* conc. of cholesterol from standard curve x volume of lipids resuspended in chloroform (week 2) Types of phospholipids present (thin layer chromatography) Announce Concepts Procedure Hazards Tips Clarification End































BSA standard curve <sup>(3)</sup>		
1 mg/mL BSA	di-water	Dye-binding reagent
0 µL	0.5 mL	4.5 mL
10 µL	0.49 mL	4.5 mL
20 µL	0.48 mL	4.5 mL
30 µL	0.47 mL	4.5 mL
40 µL	0.46 mL	4.5 mL
50 µL	0.45 mL	4.5 mL
60 µL	0.44 mL	4.5 mL
<ul> <li>Set-up your other experir volume, adjust to 500 µL</li> <li>Let tubes sit for 5 minute</li> <li>Redo any experiments th region (10 µg – 60 µg of more samples</li> </ul>	nental readings at the with di-water, add 4.4 s at room temperatur lat are outside of star BSA) → if absorbanc	e same time (protein sample 5 mL of dye-binding reagent) te then read at 595 nm indard curve interpolation te is too low, simply add









Chapter 10BC Before the lab period, you should have:
<ul> <li>✓ Completed your Pre-lab Write-up and submit on Gradescope</li> <li>✓ Title, purpose and procedures</li> <li>✓ Remember to include:</li> <li>✓ Tables for</li> <li>✓ Cholesterol determination</li> <li>✓ Dye-binding assay results</li> <li>✓ Prepare prompts for when you need to weigh and record the masses of different tubes, volumes (lyophilized samples, lipid extraction steps, and sample prep for cholesterol determination assay, volumes before lyophilization)</li> <li>At the end of lab, you should have:</li> </ul>
<ul> <li>✓ Turned in to your TFs</li> <li>✓ Dialyzed 75% Washed Supernatant in 15 mL conical tube</li> <li>✓ 25% aliquot remainder of Washed Membrane</li> <li>✓ 25% aliquot remainder of Unwashed Membrane</li> <li>✓ 25% aliquot remainder of Washed Supernatant (back-up/extra)</li> </ul>
Announce Concepts Procedure Hazards Tips Clarification End