

Announcements

- Laboratory Report Template for Chapter 10 was posted Thr 4/17 10 AM
- Chapter 10 Laboratory Report is due week of Apr 24 – 29
- If you need Acquastain and/or PAS-stained TF gels, email or Slack Dr. Ho and he will send via Slack.
- Next week Monday meets on Wednesday; Monday Discussion for Ch11 and B1 will meet at their normal times that day.

APRIL	Lecture 33	14		15	Lecture 34	16		17	Lecture 35	18
	DISCUSSION CH 10E				10D RUN THRU			DISCUSSION CH 10E	DISCUSSION CH 10E	
	CH 10D LAB			CH 10D LAB			CH 10E LAB	CH 10E LAB		
	Review Session				TOLAN EXAM 4					
		21		22	Lecture 36	23		24	Lecture 37	25
	PATRIOT'S DAY HOLIDAY				DISCUSSION CH 11C			DISCUSSION CH 11C	DISCUSSION CH 11C	
					CH 10E LAB			CH 11C LAB	CH 11C LAB	
					MONDAY SCHEDULE 11 RUN THRU			POST-LAB 11 DUE	POST-LAB 11 DUE	
								POST-LAB 10 DUE	POST-LAB 10 DUE	
	Lecture 38	28		29	Lecture 39	30		1		2
	CH 11C LAB			CH 11C LAB				LAST DAY OF CLASSES	STUDY PERIOD	
	POST-LAB 11 DUE			POST-LAB 11 DUE						
	POST-LAB 10 DUE			POST-LAB 10 DUE						
MAY	5/4-midnight	5		6		7		8		9
	Proj-3 DUE				FINALS			FINALS		FINALS
	FINALS									
	Review Session				TOLAN EXAM 5 (8-11 am)					

Announce Concepts Procedure Hazards Tips Clarification End

Chapter 10E: Lipids and Membranes

Objectives

- To probe for the presence of a glucose transporter 1 (GLUT1) in prepared fractions using a GLUT1 specific antibody (Ab)
- Assess if GLUT1 is an integral or peripheral membrane protein
- To probe for the presence of different concentrations of His-tagged KHK using a His-tag antibody

Procedures

- To complete two **immunoblot assay** using either a glucose transporter 1 (GLUT1) specific or His-tag **primary antibody**
- Visualize location of GLUT1 and his-tagged KHK with a **horseradish peroxidase conjugated secondary antibody** and **luminol with enhancers (ECL substrate)**

Announce Concepts Procedure Hazards Tips Clarification End

Membrane fractions isolated from RBCs in Ch 10

In which fraction(s) would you “expect” to observe the integral membrane protein GLUT1?

• Unwashed membrane



- **Integral and peripheral proteins**
- Membrane is intact. All membrane associated proteins are present in this fraction

• Washed membrane



- **Integral proteins**
- Peripheral proteins removed with high salt membrane wash buffer

• Washed supernatant



- **Peripheral proteins**
- Removed from membrane using high salt membrane wash buffer

Announce Concepts Procedure Hazards Tips Clarification End

Immunoblotting

- *Western Blotting – Used to identify a particular protein in a complex mixture*

- ✓ • **1st Step:** Separate proteins by SDS-PAGE
- ✓ • **2nd Step:** Transfer proteins to a membrane in same position as SDS-PAGE separation
- ☑ • **3rd Step:** Incubate membrane with blocking buffer to prevent the non-specific binding of the antibodies (TFs will do it for you)
- **4th Step:** Incubate membrane, with attached proteins, to an antibody specific for a particular protein or proteins of interest
 - Often called the **primary antibody**
- **5th Step:** Detect binding of the primary antibody using a colored or visual reaction
 - Usually requires a **secondary antibody**

Announce Concepts Procedure Hazards Tips Clarification End

Immunoblotting: detection of a target protein

Recap of what you did...

- **1st Step:** Separate proteins by SDS-PAGE
- **2nd Step:** Transfer proteins from gel to a PVDF membrane





Image: www.sigmaaldrich.com

Announce Concepts Procedure Hazards Tips Clarification End

Immunoblotting: detection of a target protein

- **3rd Step:** Blocking of membrane by incubation with blocking buffer
- Prevent or reduce nonspecific binding of the primary and secondary antibodies to proteins on the membrane or to the membrane itself
- Reduce background and improve the signal-to-noise ratio




Typically, non-fat dry milk, gelatin, or BSA are frequently used as blocking agents

Announce Concepts Procedure Hazards Tips Clarification End

Immunoblotting: detection of a target protein


- **4th Step:** Probe for the protein of interest using a primary antibody (1° Ab). The major serum antibodies are made as immunoglobulins (Ig) class-G, or IgGs.
 - For this lab, will use an antibody raised against **GLUT1*** and a **His-tag antibody**



Announce Concepts Procedure Hazards Tips Clarification End


Immunoblotting: detection of a target protein

- **5th Step:** Detect 1° Ab using a 2° Ab to visualize where target protein appeared in the original gel



- Immunoglobulin recognizes and binds to 1° Ab
- Amplify signal of 1° Ab
- Allows visualization of protein of interest through reporter molecule

IB after development*:




Images: www.bio-rad.com
www.sigmaaldrich.com

Note: This is an example only!

Announce Concepts Procedure Hazards Tips Clarification End


Immunoblotting: detection of a target protein

- **5th Step (continued):** visualization of protein of interest (chemiluminescence method)



Light signal can be captured on X-ray film or by a digital image

- However, luminol emits light only weakly, so enhancers are added to the reaction to increase the signal.
- Enhancement of a luminol-based signal is commonly referred to as enhanced chemiluminescence (ECL)



IB in imager



Announce

Concepts

Procedure


Hazards

Tips

Clarification

End

How are these reagents created?



Announce

Concepts

Procedure


Hazards

Tips

Clarification

End


How are these reagents created?



Procedure: Chapter 10E

• Probing of Immunoblots:

- You will find both of your membranes pre-blocked in blocking buffer (5% milk in 1X Tris-Buffered Saline (TBS*)) in two separate containers**. Discard blocking buffer.



- Begin with the **primary antibody incubation** step.
- Use these nice plastic trays** for incubations

*50 mM NaCl, 20 mM Tris-HCl, pH 7.5

Announce Concepts Procedure Hazards Tips Clarification End

Procedure: Chapter 10E

• Probing of Immunoblots:

- Apply 15 mL primary antibody solution. We are using anti-GLUT1 antibodies made in rabbit (**Rb** α -**GLUT1**) for your membrane fractions immunoblot and 6x-His Tag Recombinant antibodies made in rabbit (**Rb** α -**His**) for the KHK immunoblot. The 1° Abs are in 5% milk/1X **TBS**.
- Incubate membranes for 50 – 55 min at room temperature on shaker.
- ***Decant your primary antibodies (separately!!) back into Falcon tubes and give the primary antibody solution back to your TFs after using it!***
- Wash membrane with 15 mL 1X Tween-20 Tris-Buffered Saline (**TTBS**)* for 5 min on shaker (*repeat two more times*)

*0.05% (v/v) Tween-20, 50 mM NaCl, 20 mM Tris-HCl, pH 7.5

Announce > Concepts > **Procedure** > Hazards > Tips > Clarification > End

Procedure: Chapter 10E

• Staining of Immunoblots:

- Incubate membranes with ~15 mL **secondary antibody solution** (2° Ab) for 50 – 55 min on shaker at room temperature. We are using anti-rabbit IgG antibodies made in goats with the enzyme, horseradish peroxidase attached covalently (**HRP** **Goat-IgG** α -**Rb-IgG** 2° Ab). The 2° Abs are in 5% milk/1X **TBS**.
- Wash membrane twice with 15 mL 1X **TTBS** for 5 min with shaking
- Wash membrane twice with 15 mL 1X **TBS** for 5 min with shaking to remove residual detergent

Announce > Concepts > **Procedure** > Hazards > Tips > Clarification > End

Procedure: Chapter 10E

• Visualization of Immunoblots:

- You will be using the Pierce™ ECL Western Blotting Substrate.
- TFs will prepare ~12 mL substrate working solution by mixing 6 mL Detection Reagent 1 and 6 mL Detection Reagent 2*.
- **Follow your TF up to the 3rd floor! One group at a time as you are ready.*
- *On the 3rd floor, you will incubate your blot with the substrate working solution for 2-5 minutes when you are at the chemiluminescent imager.*
- Remove the immunoblot from the working solution and place it on the transilluminator. Use a Kimwipe to remove excess liquid. (You may apply a clear plastic wrap on it if desired)
- Visualize blot using the ChemiDoc™ XRS+

*This working solution is stable and can be re-used to visualize multiple blots within the hour.
(*So, this 12 mL can be shared and re-used by the entire class, ask TFs*)

Announce > Concepts > **Procedure** > Hazards > Tips > Clarification > End

Hazards

- Separate waste bottle for liquid western blot solutions

Announce > Concepts > Procedure > **Hazards** > Tips > Clarification > End

CH10E Tips

- Primary antibodies is expensive! Please do not throw it away. RETURN TO YOUR TFS.
- Save the plastic trays after your western blot so that the next section can use it.
- The ECL substrate to visualize your blot will be shared with all the groups.

Announce > Concepts > Procedure > Hazards > **Tips** > Clarification > End

Chapter 10E

- Use of ECL substrate; we will NOT use the chloronaphthol described on pages 352-353 (parts E.g and E.h). Instead, see slide 21.
- Use of His6 primary Ab in addition to the GLUT1 Ab

Announce > Concepts > Procedure > Hazards > Tips > **Clarification** > End

Chapter 10E

Before the lab period, you should have:

- ✓ Completed your Pre-lab Write-up and submit on Gradescope
 - ✓ Title, purpose and procedures
 - ✓ Staining reagent recipe
 - ✓ Flow chart of process

At the end of lab, you should have:

- ✓ Imaged both of your immunoblots (upload as part of post-lab)

Announce > Concepts > Procedure > Hazards > Tips > Clarification > End

Questions?

Discussion Quiz