

# <section-header> Acnocacce and the set of your lab (week starting Thursday Feb. 26, 20 - Tuesday Feb. 26). Introduction to Creativity Project: Goals (today) Groups (announced by Friday midnight Feb 7). Due dates (Part 1 - Feb 19; Part 2 - March 19; Part 3 - May 4).

# Chapter 6: Isolation of Plasmid DNA <u>Purpose for week 1:</u>

A) Isolate two plasmid DNAs from E. coli

B) Determine DNA concentration by two methods

→ UV absorbance (Cary-60 UV/Vis)

 $\rightarrow$  Gel electrophoresis (agarose gels)

## **Chapter 8: Transcription/translation**

Procedure Hazards Tips

Determine between plasmids A & B which is pGEM3-Rel and which is pGEM4-Rel

- pGEM3-Rel (reverse orientation of ORF)
- pGEM4-Rel (ORF in correct orientation)

Concepts







### Characterization of Plasmid Prep: Determining the DNA Concentration **UV Absorbance** • Dilute 2 $\mu$ L of plasmid mini-prep with 998 $\mu$ l of TE buffer (1:500 dilution) Record $A_{260 nm}$ values from Cary-60 specs Adjust concentration as necessary so that $0.01 < A_{\rm 260} < 1.0$ These A<sub>260 nm</sub> values are by definition readings for concentration with units of O.D./mL These O.D./mL values are for what was in the cuvette Back-calculate to the concentration of the prep by multiplying by the dilution factor (i.e., multiply by 500 for this example). Example, reading of 0.04 of 1:500 dilution is 20 OD/mL. Gel electrophoresis Take another 2-7 µL (~0.1 O.D.) of plasmid mini-prep and add electrophoresis sample buffer Load onto agarose gel and separate by size. Stain for nucleic acids and compare to standards of known amounts in ng of DNA NOTE: These are different units (OD/mL and ng/µL); will need to interconvert them to compare Concepts Procedure Hazards Tips Clarification















# **Isolation of Nucleic Acids**

- Plasmid Mini-Prep
  - Add 0.4 mL of 1:1 phenol:chloroform (v/v) to your samples
    - Pull from the bottom layer of the stock bottle; it is saturated with aqueous buffer.
    - Phenol is toxic, can cause severe burns and throat irritation
    - Chloroform is an anesthetic, but can be toxic; damage to organs
    - This step MUST be done in the hood!
  - . Tightly CAP tubes!!
  - · Vortex/shake your samples vigorously for 30 sec
  - Centrifuge to separate aqueous and organic phase
  - Transfer top aqueous phase to clean labeled tube

Concepts Procedure Hazards Tips

• Discard bottom layer and all phenol:chloroform waste directly in the hood. See Instructors.

# **Precipitation of Nucleic Acids**

- Plasmid Mini-Prep (continued)
  - Add *cold* 100% ethanol to the separated aqueous fraction, mix well
  - Centrifuge for 15 minutes at 17,000 x g
  - Remove supernatant, wash pellet with cold 70% ethanol
  - Centrifuge 1-2 minutes
  - Remove supernatant and air dry be careful not to remove pellet at the same time
  - Add 35  $\mu$ L Nuclease-free Water, vortex/pipet to dissolve pellet
  - This is your final sample = mini-prepped plasmid DNA

nnounce 🧼 Concepts 🧼 Procedure 🔪 Hazards 🔪 Tips 📄



Electrophoresis Samples				
• Sample Preparation: For Each Plasmid A & B, set up the following				
	Reagent	Volun	nes	
	Plasmid DNA: 0.1 – 0.2 OD units*	2-7 µ	μL	
	6X Sample Buffer	3 μl	μL	
	Water	15 μL – (Plasm	nid DNA μL)	
	Total Volume	18 µ	ιL	
<ul> <li>You will generate a total of 2 samples (A &amp; B) in your group</li> <li>Load Gel:</li> </ul>				
<ul> <li>Run gel with another group: 2 samples + standards/group x 2 groups = 6- 8 samples per gel</li> </ul>				
<ul> <li>Standards for each gel will be pre-aliquoted with sample buffer; load entire volume aliquoted.</li> </ul>				
<ul> <li>Supercoiled DNA Marker and/or</li> <li>DNA Mass ladder (Minnesota Molecular)</li> </ul>			*calculate the your OD/mL m in the spec.	volume from easurements
Announce	Concepts Procedure Hazard	s Tips	>>> Clarificati	on End









# Hazards

Phenol/chloroform: work in hood, PPE, care to cap tubes, watch if any liquid on gloves, use 95% ethanol to wash any off

Electrophoresis: don't stick fingers in chamber, power off when handling electrodes, power off is spilled





