

## Announcements

- Chapter 6/8 Laboratory Reports due on week of Mar 20-25.
- Chapter 8C gel image will be uploaded on Slack once we finish developing them.

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## ***Chapter 9A: Protein Modification***

### Objectives

- To learn about the enzyme ketohexokinase (KHK) and its role in fructose metabolism
- To modify sulfhydryl group (thiol group) of cysteines residues in a protein using DTNB (Ellman's Reagent)
- To measure KHK activity using a "coupled assay"
- To learn how modification of specific cysteine residues affects KHK activity

### Procedures

- To use **DTNB** to *modify cysteines groups in KHK*
- To use a **coupled assay for KHK activity** to *investigate effect of DTNB on enzyme activity*
- To **vary concentration of coupling enzymes (Pyruvate Kinase)** to *investigate KHK coupled assay*

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## Various methods to study enzymes

- Testing optimal buffer conditions
  - **Chapter 3 & 4 - LDH in Tris vs KPO<sub>4</sub> buffer**
- Steady-state kinetics with various inhibitors; measure  $K_m$ ,  $k_{cat}$ ,  $K_i$ 
  - **Chapter 4**
- Point mutation studies (*later in Biochem II*)
  - **Chapter 11C**
- **Protein modifications (Chapter 9)**

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## Protein Modifications

- **USE: Probe role of particular amino acids in protein structure and function**
- Commonly used reagents are specific for one amino acid or functional group

Reagent	Residue	Detection
2-hydroxy-5-nitrobenzyl bromide (Koshland's reagent)	Tryptophan	410 nm
N-bromosuccinimide	Tryptophan	260/280 nm
Phenylisothiocyanate (Edman's Reagent)	Amino-terminal	Release of a PTH-amino
Iodoacetic acid	Cysteine	Carboxymethyl derivati
*N-ethylmaleimide (NEM)	Cysteine	Derivatives of NEM
*5,5'-Dithiobis(2-nitrobenzoic acid) (Ellman's Reagent - DTNB)	Cysteine	412 nm
Diethylpyrocarbonate (DEPC)	Histidine	240 nm
Imidates	Lysine	Derivatives of imidates
2,4,6-trinitrobenzenesulfonic acid (TNBS)	Lysine	420 nm

\*Reagents used in this experiment; although the NEM will be attached to fluorescein such that can be detected by fluorescence.

- **ADVANTAGES:**
  - Cause an easily observed chemical change in the protein
  - Can monitor amount of incorporation into protein and change in enzyme activity
- **DISADVANTAGES:**
  - Can only observe accessible residues
  - Can cause side reactions with other residues

Announce

Concepts

Procedure


Hazards

Tips

Clarification


End

## Fructose



β-D-Fructofuranose

### High Fructose Foods



- Monosaccharide
- Food additive in many processed food

<https://fodmapfriendly.com/blogpost/how-to-enjoy-fruit-with-fructose-intolerance/>

Announce Concepts Procedure Hazards Tips Clarification End

## Fructose Metabolism





Figure 16.14  
Biochemistry, Seventh Edition  
© 2012 W. H. Freeman and Company

<http://www.bu.edu/aldolase/HFI/hfiinfo/wt.html>

Announce Concepts Procedure Hazards Tips Clarification End

## Inborn Errors in Fructose Metabolism



### Clinical Correlations

#### Fructose Metabolism: Malabsorption, Essential Fructosuria & Hereditary Fructose Intolerance

- Three maladies arise due to fructose metabolism.
- Only the last one, deficiency of aldolase B, can lead to death. Even small amounts of fructose lead to severe hypoglycemia, jaundice, hemorrhage, hepatomegaly, uricemia, and kidney damage.
- Fructose can account for 30-60% of carbohydrate in the Western diet. Humans have a limited capacity to handle this sugar, and show the same ATP depletion after a sugar bolus. Use of fructose, sorbitol, and xylitol for parenteral nutrition (i.v.) has been discontinued to to liver damage.



Giltzmann et al., 1995



Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## Kethexokinase (KHK)

- KHK exists in two isoforms – KHK-A and KHK-C



### KHK-A

- Ubiquitous expression
- $K_M = 7 \text{ mM}$  (Fructose)
- Role not clearly defined

### KHK-C

- Liver, Kidneys, Intestines
- $K_M = 0.8 \text{ mM}$  (Fructose)
- Physiological role in fructose metabolism

- Exist as a dimer
- **There are many Cys residues, the role of which are unknown.**



Human KHK Structure

Trinh et. al. (2009)

Announce

Concepts

Procedure

Hazards

Tips


Clarification

End



## 5,5'-Dithiobis(2-nitrobenzoic acid)

- DTNB reacts with Cys thiol group



- Can quantitatively determine incorporation of DTNB by  $A_{412 \text{ nm}}$

Announce

Concepts

Procedure

Hazards

Tips


Clarification

End

## 5,5'-Dithiobis(2-nitrobenzoic acid)

- Can observe effect of DTNB on activity
- Can denature enzyme and compare native and unfolded states
  - Are the number of Cys residues the same?
  - If not, why?

❖ **Buried residues may be inaccessible in native state**



Announce

Concepts

Procedure

Hazards

Tips

Clarification


End

## Measuring KHK activity



FIGURE 9-9

Recall Chapter 3...



Announce

Concepts

Procedure

Hazards


Tips

Clarification

End

## The Coupled Reaction for KHK activity

- Coupling multiple enzymes in sequential reactions
- PK and LDH must be in excess of KHK so that KHK activity is rate limiting *Why?*



**Measure KHK activity rather than PK or LDH activity!**  
 Initial rate of reaction will be proportional to [KHK]

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## General overview of Chapter 9A

### Part I: DTNB Modification Reactions

- Modify exposed Cys residues on KHK with DTNB
- Use denaturants (SDS) to test for buried Cys
- Use substrates to test for protection of potentially other important Cys
  - *Is there a difference in DTNB readings among these reagents? Why?*

### Part II: KHK Activity Assay

- Test the activity of four samples of DTNB modified KHK
  - *Is there a difference in KHK activity compared to controls? Why?*

### Part III: Coupled Assay

- Test different amounts of **PYRUVATE KINASE** to check if the coupling enzyme of the coupled assay are in excess.
  - *Why is this important?*

Announce

Concepts

Procedure

Hazards

Tips

Clarification


End

## In-Class Activity: Experimental Design

## Chapter 9A: Procedure

### Part I: DTNB Modification Reactions

- 6 reactions to set up per group (triplicates) in a [96-well clear bottom plate](#)



We only have 2 plate readers so 4 groups need to share a 96-well plate

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## Chapter 9A: Procedure

- Prepare the 6 reactions according to the table below in a 1.5 mL microcentrifuge tube
- Prepare volumes needed for **~3.2\*** reactions (slightly >3) to account for pipetting errors  
 Fold excess over reaction's final conc. **TF will add!**

Volumes used in each well***	BUFFER	KHK	Cyt NAc	ATP	SDS	H <sub>2</sub> O	Volume	DTNB	Final Volume
Fold dilution from stock	2.5x	<b>10X</b>	10x	10x	10x	–	–	–	–
[Stock]	–	<b>100 <math>\mu</math>M</b>	50 $\mu$ M	10 mM	1 %	–	–	500 $\mu$ M	–
[Final] in Assay	–	10 $\mu$ M	5 $\mu$ M	1 mM	0.1%	–	–	100 $\mu$ M	–
A	100	–	–	–	–	10	200	50	250
B	100	–	25	–	–	7	200	50	250
C	100	25	–	–	–	140	250	–	250
D	100	25	–	–	–	30	200	50	250
E	100	25	–	25	–	65	200	50	250
F	100	25	–	–	25	65	200	50	250
Order of addition:	<b>2</b>	<b>6</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>1</b>			

\*Buffer is 62.5 mM Tris, pH 7.5, and 1.25 mM EDTA. This is the 2.5x stock which has concentrations in the 1x assay of 25 mM and 0.5 mM, respectively. \*\*the stock concentration of 250  $\mu$ M is approximately 8 mg/mL. \*\*\* volumes in the table are all in  $\mu$ L

**Do not add DTNB into reaction!**

\*This is due to the limited amount of KHK we have

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## Chapter 9A: Procedure

- Record the columns that your TF assigned you to load your samples on the 96 well plate. [For example, Row A-F Column 1-3](#)
- Pipet the appropriate volumes into each well. [Each well should contain 200  \$\mu\$ L except for row C, which will not have DTNB. Row C samples should have 250  \$\mu\$ L instead](#)
- Inform your TF once you have finished loading so that the next group can continue loading
- Once 4 groups have finished loading, your TF will bring your 96-well plate to the Molecular Devices plate reader located either in SCI-156 or SCI-472 [You are welcome to follow your TFs to see how the plate readers are operated](#)

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## Chapter 9A: Procedure



Press on the “drawer” button on the plate reader to open the sample tray.



Remove the plate lid and load the plate into the drawer. **Please ensure that the plate is properly aligned in the drawer before you close it!**

Announce

Concepts

Procedure

Hazards


Tips

Clarification

End

## Chapter 9A: Procedure

- Open the **SoftMax Pro** software on the computer if not already open.
- The plate readers will be set in the **ABS (absorbance) read mode** and **Kinetic read type** for a total run time of 1 h, at 30s interval, collecting data at 412 nm



- Your TF will use a multi-channel pipet to add 50  $\mu$ L 0.5 mM DTNB to ALL samples for all 4 groups (except Row C). Check that there are no bubbles before clicking on the **READ** button
- Once the run has completed, remember to export the file **by columns** and in .xls format. You can either email the data to yourself or save it on a flash drive

Announce > Concepts > Procedure > Hazards > Tips > Clarification > End

## Chapter 9A: Procedure

### Part II: KHK Activity Assay (effect of DTNB on enzyme)

- Before measuring activity of your samples in Part II, you need to quench your DTNB reaction samples.
  - Add 100  $\mu$ L of **any one** of the triplicate KHK-containing reactions in **row C, D, E and F** into tubes containing 100  $\mu$ L 100 mM  $\beta$ -mercaptoethanol. **Put tubes on ice!**
- Note, this is a 2-fold dilution of 100/250 of your KHK added to the Part I assay.
- Cover the entire 96-well plate with the plate lid, parafilm. We will keep the entire plate at 4  $^{\circ}$ C until we can measure protein concentration

Announce > Concepts > Procedure > Hazards > Tips > Clarification > End

## Chapter 9A: Procedure

### Part II: KHK Activity Assay

- Make a "2x KHK cocktail" for at least 16 200  $\mu\text{L}$ -assays:

Reagent	Volume/2x assay cocktail ( $\mu\text{L}$ )	Final Concentration
1 TEA-HCl, pH 7.4, 1 M (30x)	133	33 mM
2 KCl, 4 M (40x)	100	100 mM
3 NADH, 30 mM (100x)	40	300 $\mu\text{M}$
4 $\text{MgCl}_2$ , 1.2 M (200x)	20	6 mM
5 Phosphoenolpyruvate (PEP), 0.133 M (100x)	20	1.33 mM
6 ATP ( $K_m = 0.13 \text{ mM}$ ), 65 mM (50x)	80	1.3 mM
7 Fructose ( $K_m = 0.25 \text{ mM}$ ), 500 mM (50x)	80	10 mM
8 Pyruvate kinase, 1000 U/mL (1000x)	4	1 U/mL
9 Lactate dehydrogenase, 1000 U/mL (1000x)	4	1 U/mL
10 $\text{H}_2\text{O}$	1519	
Total volume of KHK 2xcocktail		2000

Combine these components as part of a "cocktail." Keep cold until ~5 min prior to assay.

- The cocktail is prepared at 2x in all concentrations of reagents that are held constant. This will allow you to eventually add 100  $\mu\text{L}$  of this 2X KHK cocktail to 100  $\mu\text{L}$  of appropriately diluted KHK to give 1x concentrations of reagents.
- You can prepare the 2x cocktail while waiting for your DTNB assay, but keep it on ice.
- ~5 minutes before you start measuring activity, remove from ice and bring it to room temperature


Announce > Concepts > Procedure > Hazards > Tips > Clarification > End

## Chapter 9A: Procedure

### Making serial dilutions of KHK

You have four samples of KHK from Part I

- Take a new 96-well clear bottom plate (per group)



- Prepare dilutions of each of your four quenched KHK samples
  - Pipet 192  $\mu\text{L}$  10 mM HEPES pH 7.4, 0.5 mM DTT (EDB) into Row A, columns 1-4 (50-fold dilution) and 100  $\mu\text{L}$  EDB into Rows B, C, and D, columns 1-4
  - Add 8  $\mu\text{L}$  of quenched KHK into A1 to A4. Mix well
  - Pipet 100  $\mu\text{L}$  from row A and aliquot into row B. Mix well
  - Do the same from B to C, and C to D.
  - Discard 100  $\mu\text{L}$  from row D

Announce > Concepts > Procedure > Hazards > Tips > Clarification > End

## Chapter 9A: Procedure



FIGURE 9-9

### Part III: Coupled Assay

- Varying [Pyruvate Kinase] now instead of [KHK]
- Make **2x PK cocktail** for **at least 12 200  $\mu\text{L}$ -assays**:

	Reagent	Volume ( $\mu\text{L}$ )/assay	Final Concentration
1	TEA-HCl, pH 7.4, 1 M (30x)	100	33 mM
2	KCl, 4 M (40x)	75	100 mM
3	NADH, 30 mM (100x)	30	300 $\mu\text{M}$
4	MgCl <sub>2</sub> , 1.2 M (200x)	15	6 mM
5	Phosphoenolpyruvate (PEP), 0.133 M (100x)	15	1.33 mM
6	ATP ( $K_m = 0.13 \text{ mM}$ ), 65 mM (50x)	60	1.3 mM
7	Fructose ( $K_m = 0.25 \text{ mM}$ ), 500 mM (50x)	60	10 mM
8	Lactate dehydrogenase, 1000 U/mL (1000x)	3	1 U/mL
9	Ketohehexokinase (KHK), 5 $\mu\text{M}$ (50x; 1.5 U/mL)(from rxn C)	60	0.03 U/mL
10	H <sub>2</sub> O	1082	
	Total volume of <b>PK 2xcocktail</b>	1500	

Combine these components as part of a "cocktail."  
Keep cold until ~5 min prior to assay.

- You can prepare the **PK 2x cocktail** while waiting for your DTNB assay, but keep it on ice
- ~5 minutes before you start measuring activity, remove from ice and bring it to room temperature

Announce

Concepts

Procedure

Hazards

Tips


Clarification

End

## Chapter 9A: Procedure

### Making serial dilutions of PK

- Use the same plate for Part II, use rows A to D, column 5 to 7



- Prepare dilution of Pyruvate Kinase
  - Pipet 4  $\mu\text{L}$  of stock PK and add to 1 mL EDB
  - Add 133  $\mu\text{L}$  of diluted PK into Row A columns 5-7
- Pipet 100  $\mu\text{L}$  of EDB into wells in Rows B to D, columns 5-7
- Remove 33  $\mu\text{L}$  of samples from Row A and add to Row B. Mix well
- Do the same from B to C, and C to D.
  - Discard 33  $\mu\text{L}$  from row D

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End



## Chapter 9A: Procedure

### Part II: KHK Activity Assay & Part III Effect of lower PK on Coupled Assay

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Cyan: KHK reaction product assay  
Magenta: PK dilution assay

**Yellow part:** Aliquot 230  $\mu\text{L}$  of the 2xcocktail means for the KHK assay into each well in the yellow portion. This will allow use of the multichannel pipette to transfer liquid from each set of yellow wells twice, into two different sets of wells in the cyan area. In other words, in wells **B5-H1** you can remove 100  $\mu\text{L}$  to wells **A1-D1**, and again 100  $\mu\text{L}$  into **A2-D2**. Then **E2-H2** can be used to start reactions in **A3-D3** and **A4-D4**.

**Red part:** Similar to the yellow part, aliquot 320  $\mu\text{L}$  of the PK 2xcocktail into each well. When starting the reaction, they will use the multichannel pipette to transfer 100  $\mu\text{L}$  liquid from the bottom **E-F** wells to the top **A-D** wells in the magenta area.

- Put 2X cocktails in reservoir wells in Rows E–H (230 & 320  $\mu\text{L}$  in each well).
- Bring your 96-well plate and multi-channel pipette and tips to the plate reader
- The plate readers will be set in the **ABS (absorbance) read mode** and **Kinetic read type** for a total run time of 5-10 min, 10 s interval, 340 nm for data collection
- Once you are ready to begin, using 4 tips on the multi-channel pipet, aliquot 100  $\mu\text{L}$  of 2x cocktail into each well (A–D) of each column, and click on the **READ button**
- Again, Once the run has completed, remember to export the file **by columns** and in .xls format. You can either email the data to yourself or save it on a flash drive

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## CH 9A Hazardous chemical waste

- Leftover buffers, enzymes, and all DTNB assays go into the mainstream waste (large carboy in the fume hood)

Announce

Concepts

Procedure

Hazards

Tips

Clarification

End

## CH 9A Lab Tips

- Work together in your preparation and setting up of the DTNB reactions in the 1.5 mL microcentrifuge tubes. Split up the work if you need to. (Important to always check the pipet tips volume match the pipet you are using!)
- Be efficient in the loading of wells on the 96-well plates since 4 groups are sharing a plate
- Check for bubbles. If necessary, use a needle to remove it.
- Prepare your reagent for quenching in advance in 1.5 mL microcentrifuge tubes. Add KHK samples to these pre-assembled tubes
- While waiting for DTNB reactions to finish, start preparing cocktail for activity assay and coupled assay – keep on ice until ~5 min before starting assay, leave at room temperature
- Remember to save your samples and your data!

Announce > Concepts > Procedure > Hazards > **Tips** > Clarification > End

## Chapter 9A Clarification

- The KHK prep is less concentrated than suggested in the lab manual. It is 100  $\mu\text{M}$  instead of 250  $\mu\text{M}$ , so it is a 10X stock. You will need to re-calculate your recipes by adding more KHK and less water. Also, when making up your recipe, use 3.2-times, instead of 3.5 times what is required for 10  $\mu\text{M}$ . To make is simpler and assure you have enough KHK for all reactions, just add 75  $\mu\text{L}$  of your 10X KHK sample to each of the 4 reactions with KHK (C, D, E, &F).
- The 2x KHK Cocktail will use a mixture of pyruvate kinase and LDH. Use 8  $\mu\text{L}$  so you don't have to entirely re-calculate your recipe.
- The serial dilution of Pyruvate Kinase will be a 4-fold series. See slide 25. You will make a 500-fold dilution (4  $\mu\text{L}$  into 1 mL EDB) and use 133  $\mu\text{L}$  in each of 3 wells of row A, then do the serial dilutions.

Announce > Concepts > Procedure > Hazards > Tips > **Clarification** > End

## Chapter 9 Week 1

### Before the lab period, you should have:

- Completed your Pre-lab Write-up and submit on Gradescope
  - Title, purpose and procedures
- Remember to include:
  - Volumes of stock reagents in each DTNB reaction tube
  - Table for DTNB reactions ([Part I](#))
  - Recipes for KHK activity assay ([Part II](#)) and Coupled assay ([Part III](#))

### At the end of lab, you should have:

- Completed DTNB modification reactions (Part I) downloaded from spec.
- Completed KHK activity assays (Part II) downloaded from spec.
- Completed Coupled Assays (Part III) downloaded from spec.
- Saved and turned in to TFs the 96-well plate from DTNB reactions and any remaining quenched aliquots for Week 2
- Record columns that your TF assigned you to load your samples on the 96 well plate for Part I
- Saved all your data

[Announce](#)[Concepts](#)[Procedure](#)[Hazards](#)[Tips](#)[Clarification](#)[End](#)

## Questions?

## Discussion Quiz