


FV10i software

Menu with detail walkthrough

FV10i software startup


- Press and release power button on front of the 10i
 - BEFORE turning on: Left button **red**, right button off
 - Warming Up: Left button **green blink** , right button **solid red**
 - Ready to start: Left button **green steady**
- Log in to the “Confocal User” account, there is No PWD
- Double click the icon on the desktop to launch the software.
 - The main unit warms up for the first three minutes after being turned on. Do not activate the software during this period.
- Enter [User ID] and [Password], and then click the <OK> button.
- This software launches and the two Fluoview windows appear

Tool Bar (list of Tabs)




1. Sample Setting
Fluorophore selection
2. Acquire Map Image
Overview of sample
3. Observe
Actual scanning
4. Review
Examine in FV explorer
5. Exit
Quit the program
- 0-Property-
ADVANCED USE ONLY
These settings effect all of the users
on your lab account!

Observe



- A. Image Acq Mode Setting
- B. Live View display
- C. Setting various parameters of scan
- D. Channel sensitivity and laser intensity
- E. Registering image acq points
- F. Starting the acq.

Observe A



1. Image Acq Mode
2. Select Bin
3. Navigation- step by step walkthrough
4. Show or hide map image channels
4. Map Image of active channels
4. Zoom (of map image)

Observe B



Repeat-start scanning (Peach when active)

Individual channels (toggle on and off)

Subscan, split channels, LUT adjust
Zoom changes mag with slider (lens with button)

Rotation

Information at the bottom of the image windows with Image size, pixel resolution, X and Y coordinates and intensity

Observe C



Theoretical Information:\size of pixels in X, Y and Z (thickness of optical slice)

Optimize

zoom: Nyquist for X/Y


Z stack: Nyquist for Z

Depth- Z location of the stage, mark points and step size for Z scanning

Other Parameter

- image size
- Averaging
- Confocal Aperture (pinhole) size in airy units
- Cross Talk Correction: sequential scanning to minimize bleed through of signals

Observe D



Sensitivity-For a fluorescence channel, the detector sensitivity (voltage applied to the photomultiplier) is adjusted. When at 50%, it corresponds to 625V. For a phase contrast channel, the target intensity value is adjusted. When this field is set to 50%, it corresponds to a target intensity value of 2048.

Laser

Percent of maximum laser power per laser

When more than two channels are active the laser must be highlighted to redraw that channel in focus mode.

Observe F



Z: [Number of slices for which images have already been acquired] / [Number of all slices in Z Stack]

A: [Number of areas for which images have already been acquired] / [Number of all areas in Multi Area]

R: [Number of repeats of image acquisition already completed] / [Total number of repeats for Time Lapse]

Save Image Folder: autosave path for images (should be your lab and subfolder)

File name: autoincremented base file name will have 01 added to second capture


One shot: captures a single z plane, timepoint area with whatever channels are active

Data File Management



Copy Data from your user folder in WINDOWS explore by selecting the .oif and .oif.files folder

(eg notrash was the name of the acquisition. Notrash.oif is a small 14kB file that will launch the full resolution version of the files stored in notrash.oif.files (folder exploded at right)

If you copy only the .oif file from Fluoview explorer, you will have a low res thumbnail only



notrash3



notrash3

(Fiji Is Just) ImageJ 2.0.0-rc-68/1.52i; Java
1.8.0_66 [64-bit]; Windows 10 10.0; 364MB
of 18178MB (2%)

Title: notrash_03.oif

Width: 212.1333 microns (1024)

Height: 212.1333 microns (1024)

Size: 44MB

Resolution: 4.8272 pixels per micron

Voxel size: 0.2072x0.2072x0.799 micron³

ID: -5

Bits per pixel: 16 (unsigned)

Display range: 307 - 4079

Frame: 8/22 (c:1/3 z:8/22 - Series 1)

Frame interval: 1 sec

No threshold

Magnification: 0.50

ScaleToFit: false

Uncalibrated

Color calibration: 12 11 1 (1526 264)

